

Product data sheet

www.dinolite-uk.com

Dino-Lite MS33W stand

Round metal base with flexible arm and holder for the Dino-Lite. The MS33W has the possibility to vary the height (appr. 15 mm) with a turning knob in the base.

The MS33W desktop stand holds the Dino-Lite with a flexible gooseneck arm that can be bent and twisted into many positions. An additional feature of the MS33W is the fine adjustment knob in the base of the stand. Using this knob, you can focus the Dino-Lite in situations where it may be too difficult to move the adjustment dial on the microscope itself.



Technical data

Model	MS33W stand
Unit dimensions	flexible arm: 30cm