

Product data sheet

www.dinolite-uk.com

Dino-Lite MS09B stand

The MS09B is an ultra-portable mini stand that holds the Dino-Lite in place whenever you need to.

Two different bases are included; a horseshoe shaped base to examine an object below, and a plate style base on which an object can be placed to examine. The plate style base also has grooves to help orientate objects placed on top of it. Both bases have rubber grips to keep the stand in place and protect the desktop surface.



The MS09B is a great little accessory to enhance the Dino-Lite's microscopic portability solution.

Technical data

| Model | MS09B |
|----------------------------|----------|
| Height | 13,5 cm |
| Weight with horseshoe base | 81 gram |
| Weight with plate base | 125 gram |
| Package weight | 186 gram |