

Product data sheet

www.dinolite-uk.com

Dino-Lite AM7013MT

The Dino-Lite AM7013MT has a 5 megapixel sensor for crystal clear images, even under higher magnification. Because of the very high resolution, the AM7013MT is the ideal solution for displaying microscopic images on large screens or using the microscopic images for printing.



Technical data

Model	AM7013MT
Resolution	5 megapixel (2592x1944)
Magnification	10x – 60x & 200x
Interface	USB 2.0, High speed USB transmission
Nr. of LED's	8
LED on/off switchable	√ (by software)
LED color	White
Material	Full aluminum alloy body
Microtouch	√
Output	Picture, video and time-lapsed video
Video Frame Rate	up to 30 frames-per-second
Measurement	√ (line, radius, circle, 3-point circle, etc.)
Calibration	√
Filter/Diffuser	-
Compatibility	Windows XP, Vista, 7 and recent Mac OS
Software	DinoCapture for Windows included. Please note that a Motion-JPEG codec is required to display 5 megapixel resolution (not included)
Other	Ideal solution for displaying images on large screens or for creating print-ready images