

Dino-Lite MS36B stand

Square metal base with column and support providing several adaptable heights. The MS36B has a horizontal beam to enable additional flexibility in the horizontal plane.

The MS36B stand is designed for use with the complete line of Dino-Lite digital microscopes and offers many fine adjustment features.

This sturdy and versatile stand utilizes an adjustable sliding horizontal boom arm with up to 25cms of reach and 360 degrees rotational capability making it ideal for inspection applications where the ability to extend out over and across a specimen is required.

The MS36B provides smooth horizontal and vertical adjustment positioning. It is manufactured from a combination of resilient stainless steel and lightweight recycled aluminum making it very portable.



Technical data

Model	MS36B
Unit dimensions (WxLxH)	12,7 x 20,7 x 30,5 cm
Other	vertical pole can be extended with MS30X kit