

## Product data sheet

[www.dinolite-uk.com](http://www.dinolite-uk.com)

### Dino-Lite AM4013MZT4

The AM4013MZT4 with a high magnification of 400~470x reveals even the smallest details on all kinds of objects. The built-in, switchable polarizer makes the AM4013MZT4 the right choice for working with high glare materials such as jewelry, plastics and metals. The aluminum housing offers enhanced protection and durability and has the best look and feel.



#### Technical data

<b>Model</b>	<b>AM4013MZT4</b>
Resolution	1,3 megapixel (1280x1024 pixels)
Magnification	400~470x
Interface	USB 2.0, High speed USB transmission
Nr. of LED's	8
LED on/off switchable	√ (by software)
LED color	White
Material	Full aluminium alloy body
Microtouch	√
Output	Picture, video and time-lapsed video
Video Frame Rate	up to 30 frames-per-second, 15fps at 1.3M
Measurement	√(line, radius, circle, 3-point circle, etc.)
Calibration	√
Filter/Diffuser	Built-in polarizer
Compatibility	Windows XP, Vista, 7 and recent Mac OS
Software (Included)	Windows: DinoCapture, Mac: DinoXcope
Other	The polarizer can be switched on/off or can be adjusted to offer full or half polarization effect